



HYPOXIC WORKSTATION



01 ACCURACY

- Precise Control of O₂ and CO₂
 - ✓ Range of O₂ Control: 0.0% - 20.9%, in 0.1% increments
 - ✓ Range of CO₂ Control: 0.03% - 20.0%, in 0.1% increments
- Precise Control of Temperature and Humidity
 - ✓ Range of Temperature Control: 5°C above ambient up to 45°C, in 0.1°C increments (Optional high/low temperature control chamber)
 - ✓ Range of Humidity Control: Room humidity to 75%, in 0.1% increments
- Advanced Airway Design
 - ✓ Ensures a high level of homogeneity in the workstation's gas distribution



- Accurate Gas Control System
 - ✓ High precision O₂ and CO₂ sensor
 - ✓ Precise mixing and control of the internal atmosphere

02 CLEANLINESS

- High efficiency HEPA filtration system.
- The top, front and bottom panels are moulded from a single piece of acrylic to reduce condensation and microbial growth.
- The clear acrylic front panel of the workstation effectively avoids front panel blur caused by vapour condensation.



03 SAFETY

- A double pressure relief system consisting of a pressure relief tank and various pressure relief valves ensures the safety of the workstation.




Pressure Relief Tank

04 CONVENIENCE

- Direct-access operation system
 - ✓ The sleeves are made of a soft, skin-friendly material and an efficient vacuum pump direct access to the workstation only 8-12 seconds.
 - ✓ Can be operated directly with lab gloves.
 - ✓ Glove port caps can be "parked" in special holders and do not take up valuable working space.



- One-touch calibration for O₂ and CO₂ sensor
 - ✓ The O₂ and CO₂ sensors can be automatically calibrated to guarantee the accuracy of gas concentration detection.



- Thoughtful lighting configurations: Daylight lamps and fluorescent detection lamps



- Programmable four-stage hypoxic cycle
 - ✓ Can stimulate the sample by change of O₂ and CO₂ concentration automatically.

Hypoxic Cycle Settings						
O ₂	CO ₂	N ₂	2020.6.22		20:20:24	
PARAMETER SETTING:			STATE: WAITING			
	O ₂ (%)	CO ₂ (%)	TIME(MIN)	O ₂ VALUE		MINS
STEP1	10.0	5.0	2	11.0 %	ACT1	0
STEP2	8.0	5.0	2		ACT2	0
STEP3	5.0	5.0	2	CO ₂ VALUE	ACT3	0
STEP4	2.0	5.0	2	4.9 %	ACT4	0
CYCLES	3			CYCLES		0
			START		BACK	

- During the interlock purge the user can select the oxygen concentration required in the interlock


Hypoxic Workstation Control Panel						
O ₂	N ₂	CO ₂	2022			
O ₂ Set Point	8.0 %	O ₂ Cur Value	1%	2%	3%	4%
H ₂ Set Point	0.0 %	CO ₂ Cur Value	5%	6%	7%	8%
Humidity Set	55 %	Humidity Cur	9%	10%	Normal	Stop
Running State: Waiting						
Start		O ₂ Cal	CO ₂ Cal	Interlock	Light-Off	

- An optional built-in air pump and filter are used to provide the oxygen supply, avoiding the need for a compressed air cylinder
- USB data logging system holds 3 months of continuous data storage
- Optional anaerobic system: hypoxic and anaerobic system in one workstation: a very low O₂ level of <5 PPM O₂ can be achieved


05 CUSTOMISED PRODUCTS




Dual Hypoxic Workstation DY-D (Large Interlock)




Hypoxic Workstation For Animal & Cell Research DY-WC-D (L: Animal; R: Cell)



Hypoxic Workstation with a Larger Workspace DY-TM



Hypoxic Workstation with Three Large Chambers DY-T



Integration of Microscopes into a Hypoxic Glove Box

Tel: 8610-88693537

Email: sales@maworde.com

www.maworde-biotech.com

